



Product Catalogue

EuniceMed Health Care



Orthopedic support / Medical compression stocking / Footcare



NOT JUST A MOTION

Freedom to move unrestricted and safely at any time.

Don't stop. Go wild and do what you love, no matter what age you are.
That is your passion calling you.

Take the adventure and witness your unlimited potential thru self-challenge.

Create happy moments with your family and friends.

EuniceMed believes the true spirit of Motion is to improve life quality.




NOT JUST A MOTION



Design and Produce calculated.

EuniceMed is a brand launched by the Comfort Plus Corporation. A Taiwan-based manufacturer focused on producing medical devices such as Orthopedic Supports, Medical Compression Stockings, Medical Grade Silicone Insole and foot care accessories for more than two decades. The solid background behind EuniceMed stands for authentic quality of comfort, durability, safety, and an enjoyable wearing experience.

We believe there is always a deeper context in one's motion, such as freedom, passion, challenging spirit, happiness, etc. Therefore, we aim to keep your motion unconstrained from discomfort and live freely every moment.



Orthopedic support

We aim to keep your motion unrestricted from discomfort and live freely every moment.

Orthopedic support is a device a person wears to assist with their daily function, reduce pain with movement, or improve efficiency in a task they want to return to. It provides firm support, aligns, prevents, or enhances the function of the movable parts of the body.

Supports can be used around the knee, ankle, elbow, wrist, and back.

EuniceMed Product Category

Compression stockings

Compression stockings help relieve symptoms and slow the progress of varicose veins. They improve circulation and are a mainstay of treatment for varicose veins that are causing symptoms.



Insoles and Footcare

The insole and heel pad are designed for better comfort and fitting. Made of 100% Medical Grade Silicone, it helps to relieve heel pain and discomfort from sore hurting heels, sore feet, plantar fasciitis, bone or aching heel spurs, Achilles, Edema, treat stress fractures, and sprained feet when walking or running.

Relieving your tired throbbing feet, giving you greater Arch Support. protect and cushion your sore heels from Shock, Bruising & Impact.





Orthopedic Support

Solution for joint discomfort.



Breathable: Breathable material brings comfort.



Far infrared ray of Nano bamboo charcoal: Nano bamboo charcoal with far infrared ray helps improve circulation and keep users healthy.



Enhanced support: Additional stay to enhanced support and prevent injuries when users moving around.



Silicone/Gel pad: To decompress and lightly massage joints.



Strengthened straps provide comfortable compression: Help stabilize joints and prevent injuries.



Patent product: Innovative patent product provide customers excellent protection and all-embracing performance with palmary comfort.

Orthopedic Support



CPO-1304

Elbow Support IU+

Elbow support with silicone pads

Feature

- Gel pad provides extra compression to stabilize the injured area.
- Provides reliable support during stressful activities.
- Supports forearm tendons & muscles, golfer's or tennis elbow relief.

Application

Epicondylitis
Tennis / Golfer's elbow

Material

Nylon, Spandex, Gel pad.

Quantity

1 PC

M	20 - 25	cm	8 - 10	in
L	25 - 30	cm	10 - 12	in
XL	30 - 35	cm	12 - 14	in





CPO-1305

Elbow Support IU

Bamboo elbow Support

Feature

- Support weak elbow and strengthen the muscle.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of elbow joint.
- Ergonomically designed for a full range of motion, protection, and compression.

Application

Weak elbow joint
Elbow pain

Material

Nylon, Spandex, Bamboo charcoal fiber.

Quantity

1 PC

S	15 - 20	cm	6 - 8	in
M	20 - 25	cm	8 - 10	in
L	25 - 30	cm	10 - 12	in
XL	30 - 35	cm	12 - 14	in





CPO-1401

Wrist Support IU

Wrist Support

Feature

- Help move freely with more support and less discomfort.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of wrist joint.
- Advanced knitting technology and high-breathable materials.

Application

Weak wrist joint
Wrist pain

Material

Nylon, Spandex, Bamboo charcoal fiber.

Quantity

1 PC

S	14.5 - 16.5	cm	5 ¾ - 6 ½	in
M	16.5 - 18.5	cm	6 ½ - 7 ¼	in
L	18.5 - 21.5	cm	7 ¼ - 8 ½	in
XL	21.5 - 23.5	cm	8 ½ - 9 ¼	in





CPO-1603

Knee Support IU

Bamboo knee support

Feature

- Advanced knitting technology with high-breathable materials.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of knee joint.
- To hold the posture correctly and provide proper compression on the knee joint.

Application

Weak knee joint
Knee pain

Material

Nylon, Spandex, Bamboo charcoal fiber.

Quantity

1 PC

S	28.5 - 33.5	cm	11 ¼ - 13 ¼	in
M	33.5 - 39	cm	13 ¼ - 15 ½	in
L	39 - 44.5	cm	15 ½ - 17 ½	in
XL	44.5 - 50	cm	17 ½ - 19 ¾	in



Ankle and Shin Sleeve

Ankle and shin sleeve



CPO-1604

Feature

- Provide suitable compression on legs, and improve the performance of muscle strength.
- Provides firm support help to increase the stability of ankle joint.
- Open-toe design.

Application

Tired legs and ankles
Weak shin

Material

Nylon, Lycra.

Quantity

1 PR



S	28 - 31	cm	11 - 12 ¼	in
L	31 - 35	cm	12 ¼ - 13 ¾	in

Shin Sleeve

Shin sleeve



CPO-1605

Feature

- Provide suitable compression on legs, and improve the performance of muscle strength.
- Helps relieve discomfort from an injured or weak shin.
- Easy on and off.

Application

Weak shin
Tired leg

Material

Nylon, Lycra.

Quantity

1 PR



S	28 - 31	cm	11 - 12 ¼	in
L	31 - 35	cm	12 ¼ - 13 ¾	in

Knee Support Auxi+

Knee support with silicone pad



CPO-1606

Feature

- Anatomically shaped compression sleeve perfectly gives support and relieves aching knee joint.
- Bilateral dual flexible spiral stay stabilizers strengthen knee joint stability and reduce knee load.
- Patella tendon pad is designed perfectly to fit the kneecap curve to reduce discomfort on the patella.

Application

Osteoarthritis
Patella instability
Weak knee joint
Knee pain

Material

Nylon, Spandex, Silicone pad, Steel stay.

Quantity

1 PC



S	33 - 36	cm	13 - 14	in
M	36 - 38	cm	14 - 15	in
L	38 - 41	cm	15 - 16	in
XL	41 - 46	cm	16 - 18	in
XXL	46 - 52	cm	18 - 20 ½	in

Knee Support Auxi

Knee support with spiral stays



CPO-1609

Feature

- Maximum pressure relief.
- A breathable compression sleeve allows more stability when moving.
- Bilateral dual flexible spiral stay stabilizers strengthen the knee joint and reduce knee load.

Application

Osteoarthritis
Weak knee joint
Knee pain

Material

Nylon, Spandex, Steel stay.

Quantity

1 PC



S	33 - 36	cm	13 - 14	in
M	36 - 38	cm	14 - 15	in
L	38 - 41	cm	15 - 16	in
XL	41 - 46	cm	16 - 18	in

Orthopedic Support

Knee Support UE

Knee support with gel pad



CPO-1613

Feature

- Circular contoured pad stabilized the patella.
- High-elastic tight and breathable compression sleeve to reduce the pressure on knee.
- Anatomically shaped compression sleeve help to hold the posture correctly.

Application

Patella instability
Weak knee joint
Knee pain

Material

Nylon, Spandex, Gel pad.

Quantity

1 PC



S	28.5 - 33.5	cm	11 ¼ - 13 ¼	in
M	33.5 - 39	cm	13 ¼ - 15 ½	in
L	39 - 44.5	cm	15 ½ - 17 ½	in
XL	44.5 - 50	cm	17 ½ - 19 ¾	in

Ankle Support IU

Bamboo ankle support



CPO-1703

Feature

- Protects the weak ankle joint and relieves discomfort from injured.
- Provide proper compression on the ankle joint.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of ankle joint.

Application

Ankle pain
Weak ankle joint

Material

Nylon, Spandex, Bamboo charcoal fiber.

Quantity

1 PC



S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in

Ankle Support IU+

Ankle support with gel pads



CPO-1704

Feature

- Bilateral gel pads provide soft compression and stabilize the ankle.
- Help reduce swelling and discomfort.
- Provide suitable compression on the ankle joint.

Application

Ankle instability
Ankle pain or weakness
Ankle sprain

Material

Nylon, Spandex, gel pad.

Quantity

1 PC



S	15 - 20	cm	6 - 8	in
M	20 - 25	cm	8 - 10	in
L	25 - 30	cm	10 - 12	in
XL	30 - 35	cm	12 - 14	in

Tennis Elbow Support

Tennis elbow support



CPO-2301

Feature

- Silicone pad provides extra compression to crucial muscle points, to stabilize the injured area.
- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.

Application

Tennis elbow
Golfer's elbow

Material

Neoprene, Nylon, Silicone pad

Quantity

1 PC



One Size

Golf / Tennis Elbow Support Strap

Tennis / Golf elbow support



CPO-2302

Feature

- Exclusive innovative patent for excellent protection and comfort.
- A medical-grade silicone pad is designed to alleviate the forearm's stress and fatigue while protecting the tendon from further discomfort.
- Lightweight design helps to improve performance.

Application

Tennis elbow
Golfer's elbow

Material

Neoprene, Nylon, Silicone pad.

Quantity

1 PC

One Size



Wrist and Thumb Support

Wrist and thumb support



CPO-2401

Feature

- Mild compression support with stabilizing strap.
- Increased comfort.
- Help to relieve stiffness palm or discomfort wrist.

Application

Stiffness hand
Wrist weakness
Wrist sprain

Material

Neoprene, Nylon, Spandex.

Quantity

1 PC



S	14.5 - 16.5	cm	5 ¾ - 6 ½	in
M	16.5 - 18.5	cm	6 ½ - 7 ¼	in
L	18.5 - 21.5	cm	7 ¼ - 8 ½	in
XL	21.5 - 23.5	cm	8 ½ - 9 ¼	in

Thumb Support AI+

Thumb splint



CPO-2402

Feature

- Effective stabilization.
- Use for De Quervain's syndrome.
- Aluminum stay for thumb immobilization.

Application

De Quervain's syndrome
Thumb sprain

Material

Neoprene, Nylon, Spandex,
Aluminum stay.

Quantity

1 PC



S / M	14.5 - 18.5	cm	5 ¾ - 7 ¼	in
L / XL	18.5 - 23.5	cm	7 ¼ - 9 ¼	in

Wrist Brace

Enhanced wrist splint



CPO-2403

Feature

- For carpal tunnel syndrome use.
- Two aluminum stays help to stabilize the wrist, reduce discomfort, and limit movement to avoid re-injured.
- Easy and adjustable strap allows for each individual to customize pressure.

Application

Carpal tunnel syndrome
Wrist sprain
Colles fracture

Material

Neoprene, Nylon, Aluminum.

Quantity

1 PC



S	14.5 - 16.5	cm	5 ¾ - 6 ½	in
M	16.5 - 18.5	cm	6 ½ - 7 ¼	in
L	18.5 - 21.5	cm	7 ¼ - 8 ½	in
XL	21.5 - 23.5	cm	8 ½ - 9 ¼	in

Orthopedic Support



CPO-8401

Wrist Thumb Brace

Wrist brace with thumb support

Feature

- For carpal tunnel syndrome use.
- For De Quervain's syndrome use.
- Five aluminum stays limited the range of motion of thumb, and enhanced the stability of wrist.
- Adjustable velcro with strip helps to fit the wrist joint.

Application

Carpal tunnel syndrome
De Quervain's syndrome
Wrist sprain
Colles fracture

Material

Nylon, Polyester, Aluminum

Quantity

1 PC

Right / Left

S	14.5 - 16.5	cm	5 ¾ - 6 ½	in
M	16.5 - 18.5	cm	6 ½ - 7 ¼	in
L	18.5 - 21.5	cm	7 ¼ - 8 ½	in
XL	21.5 - 23.5	cm	8 ½ - 9 ¼	in





CPO-2404

Wrist Support AG

Wrist wrap

Feature

- Very comfortable for protection, and injury prevention.
- Enhance support for mild joint pain relief.
- Strengthens weak muscle via adjusting straps.

Application

Wrist pain
Wrist weakness

Material

Neoprene 65%, Nylon 35%

Quantity

1 PC

One Size





CPO-2405

Wrist Support AG+

Wrist wrap

Feature

- Enhances support for mild joint pain relief and protect weak muscle.
- Premium breathable neoprene provides high stability to the wrist joint.
- Easy adjusting the support freely to the most comfortable wrapping states.

Application

Wrist weakness
Wrist pain

Material

TPE, Nylon, Polyester.

Quantity

1 PC

One Size





CPO-6403

Wrist Support IU+

Wrist support

Feature

- Adjustable strap provide perfect compression and stability for wrist joint.
- Strengthen weak wrist joint.
- Protect and stabilize the injured wrist joint.

Application

Wrist pain
Wrist weakness

Material

Nylon, Rubber, Polyester.

Quantity

1 PC

One Size



Thumb Support AI

Thumb splint



CPO-6404

Feature

- Stabilized the thumb to prevent re-injured.
- Use for De Quervain's syndrome.
- Mouldable aluminum splint to limit the thumb motion and reduce discomfort.

Application

De Quervain's syndrome
Thumb sprain

Material

Cotton, Rayon, Rubber, Aluminum

Quantity

1 PC



S	13 - 17	cm	5 - 6 ¾	in
L	17 - 25	cm	6 ¾ - 9 ¾	in

Knee Support ST

Knee support



CPO-2601

Feature

- Very comfortable.
- Mild joint pain relief.
- Adjusting straps to strengthen weak muscles and knee.

Application

Knee pain
Knee weakness

Material

Neoprene, Nylon.

Quantity

1 PC



One Size				
----------	--	--	--	--

Knee Brace GT+

Hinged knee brace



CPO-2602

Feature

- Bilateral hinged brace with straps provides maximum stability and compression to enhance faster recovery.
- The front open patella with padding stabilizes the patella and relaxes tendons, avoids pulling, and prevents muscle strain.
- Provide perfect support to the weak knee.

Application

Weak or overstressed knee
Mild strain / sprain of knee
Osteoarthritis

Material

Neoprene, Nylon, Spandex, Steel stay.

Quantity

1 PC



S	33 - 37	cm	13 - 14 ½	in
M	37 - 41	cm	14 ½ - 16	in
L	41 - 45	cm	16 - 17 ¾	in
XL	45 - 49	cm	17 ¾ - 19 ¼	in

Knee Brace GT

Knee brace with circular strip



CPO-2603

Feature

- Bilateral flexible spiral stays enhance the stability of the knee joint.
- Adjustable strips help to stabilize the knee joint.
- Front open patella with an extra silicone pad enhances the stability of the kneecap.

Application

Weak or overstressed knee
Osteoarthritis
Mild strain / sprain of knee

Material

Neoprene, Nylon, Silicone pad, Steel stay.

Quantity

1 PC



S	33 - 37	cm	13 - 14 ½	in
M	37 - 41	cm	14 ½ - 16	in
L	41 - 45	cm	16 - 17 ¾	in
XL	45 - 49	cm	17 ¾ - 19 ¼	in

Orthopedic Support



CPO-2605

Knee Support AG

Knee support with open knee

Feature

- Easy and comfortable to wear.
- Lightweight and skin-friendly.
- Enhance support for mild joint pain relief and protect weak muscle.
- For relief and mobilization.

Application

Patella instability
Knee weakness
Knee pain

Material

Neoprene, Nylon, Silicone pad.

Quantity

1 PC

One Size





CPO-2606

Patellar Knee Strap

Patellar strap

Feature

- Extra compression on patellar tendon.
- Relieves discomfort from Chondromalacia Patella.

Application

Chondromalacia Patella
Patella instability

Material

Neoprene, Nylon, Spandex, EVA.

Quantity

1 PC

One Size





CPO-8601

Patellar Knee Strap

Patellar strap

Feature

- Extra compression on patellar tendon.
- Relieves discomfort from Chondromalacia Patella.
- Eco-friendly TPE Material: Mainly made of TPE, which is non-toxic, and recyclable plastics.

Application

Chondromalacia Patella
Patella instability

Material

TPE, Nylon, Spandex, EVA.

Quantity

1 PC

Environmentally Friendly



S	25 - 34	cm	10 - 13 ½	in
Onesize	27 - 38	cm	10 ¾ - 15	in
L	29 - 42	cm	11 ½ - 16 ½	in





CPO-2607

Knee Support ST+

Knee support with patella warp

Feature

- Flexible side spiral stays strengthen the knee joint.
- Breathable Neoprene provides comfortable wear.
- The circular opening around the kneecap helps reduce pressure directly over the sore patella.

Application

Patella instability
Knee weakness
Knee pain

Material

Neoprene, Nylon, Steel stay, silicone pad.

Quantity

1 PC

S	33 - 41	cm	13 - 16	in
L	41 - 49	cm	16 - 19 ¼	in



Knee Support AG+

Knee support with open patella



CPO-2608

Feature

- An anatomically contoured pad surrounds the kneecap to stabilize the patella and relieve joint discomfort.
- Perfect for sports protection, exercise, and injury recovery.
- Spiral stays fit tightly with your knee joint; it helps to strengthen side support and maintain joint stability.

Application

Patella instability
Knee weakness
Knee pain

Material

Neoprene, Nylon, Polyester, Steel stay.

Quantity

1 PC



Normal	36.5 - 45	cm	14 ¼ - 17 ¾	in
Enlarged	45 - 51	cm	17 ¾ - 20	in

Ankle Support ST

Ankle support



CPO-2701

Feature

- Ankle sprain recovery support.
- Provide ankle stability while doing activities.
- Easy and comfortable to wear.

Application

Ankle instability
Ankle pain or weakness
Ankle sprain / strain

Material

Neoprene, Nylon.

Quantity

1 PC

S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in

Ankle Support ST+

Ankle support



CPO-2702

Feature

- Ideal for sporting.
- Prevent ankle sprains or ankle dislocations.
- Helps relieve aching muscles and injured ankle.

Application

Ankle instability
Ankle pain or weakness
Ankle sprain / strain

Material

Neoprene, Nylon.

Quantity

1 PC

S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in

Ankle Support Auxi+

Ankle support with Achilles protection pad



CPO-2705

Feature

- Easy and adjustable straps to maintain the stability of the ankle joint.
- Innovative gel pad to protect the Achilles tendon.
- Ankle joint stabilizer.

Application

Achilles tendinitis
Ankle instability
Ankle sprain / strain

Material

Neoprene, Nylon, Polyester, Gel pad.

Quantity

1 PC



S	20 - 26	cm	8 - 10 ¼	in
L	26 - 32	cm	10 ¼ - 12 ½	in

Orthopedic Support

Ankle Support Auxi

Ankle support



CPO-2706

Feature

- Provide stabilization on the ankle joint.
- Helps relieve aching muscles and injured ankle.
- Lightweight and comfortable to wear.

Application

Ankle pain or weakness
Ankle instability
Ankle sprain / strain

Material

Neoprene, Nylon.

Quantity

1 PC

S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in

Foam Cervical Collar

Foam cervical collar



CPO-3101

Feature

- To stabilize the cervical spine in neutral position.
- Velcro tape for easy adjustment.

Application

Cervical soreness
Cervical sprain / strain

Material

Cotton yarn with F/I, PU

Quantity

1 PC

S	25.5 - 33	cm	10 - 13	in
M	33 - 40.5	cm	13 - 16	in
L	40.5 - 48	cm	16 - 19	in



Philadelphia Collar

Philadelphia collar



CPO-4101

Feature

- Provides the stability and limited the range of motion of cervical vertebra.
- Anterior opening for trachea.
- Hook and loop closure for easy adjustment.

Application

Cervical spondylosis
Post-operative use for
cervical care
Cervical fracture

Material

Foam, PP, Nylon.

Quantity

1 PC

S	25.5 - 33	cm	10 - 13	in
M	33 - 40.5	cm	13 - 16	in
L	40.5 - 48	cm	16 - 19	in



Cervical Collar

Cervical collar



CPO-4102

Feature

- Adjustable-height maintaining the neck in neutral position.
- Breathable holes decrease muggy feelings.

Application

Post-operative use for
cervical care
Cervical sprain / strain

Material

PE, PVC, Nylon

Quantity

1 PC

S	33 - 38	cm	13 - 15	in
M	38 - 43	cm	15 - 17	in
L	43 - 48	cm	17 - 19	in



Ankle Brace Air

Adjustable air ankle brace



CPO-4702

Feature

- Adjustable air pad design perfectly fit on ankle joint.
- The splints in both sides to stabilize injured ankle joint from inversion and eversion.

Application

Ankle sprain / strain
Acute injuries of ankle joint
Mild swollen ankle

Material

PP, PVC, Nylon, PV.

Quantity

1 PC



SIZE	EURO	UK	US	
			MEN	WOMEN
S	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

Ankle Brace Iced

Gel ankle brace



CPO-4703

Feature

- Cold compression therapy to ankle joint by iced insert pad.
- Insert pad could be took off and put in the refrigerator before use.
- The splints in both sides to stabilize injured ankle joint from inversion and eversion.

Application

Ankle sprain / strain
Acute injuries of ankle joint
Mild swollen ankle

Material

PP, Nylon, Gel, PV.

Quantity

1 PC



One Size				
----------	--	--	--	--

Shoulder Support

Shoulder support



CPO-3102

Feature

- To stabilize and support the shoulder joint.
- Adjustable hook and loop closure for easy and convenient application.
- Can be worn on left or right shoulder.

Application

Shoulder tendonitis
Shoulder pain
Shoulder instability

Material

Polyester, Nylon, PU foam.

Quantity

1 PC

S / M	20.5 - 30.5	cm	8 - 12	in
L / XL	30.5 - 40.5	cm	12 - 16	in

Arm Sling

Arm sling



CPO-7301

Feature

- Breathable and lightweight mesh material provides comfortability.
- Adjustable velcro strap to customize fit even after putting on the sling.
- Stabilized the position of the shoulder joint and arm.

Application

Forearm fracture
Shoulder dislocation
Brachial plexus injury

Material

Polyester.

Quantity

1 PC



S	23.5 - 25.5	cm	9 ¼ - 10	in
M	25.5 - 30.5	cm	10 - 12	in
L	30.5 - 35.5	cm	12 - 14	in
XL	35.5 - 40.5	cm	14 - 16	in

Orthopedic Support



CPO-6202

Back Support GT+

Lumbar support

- 9 inches in height with four steel stays provide rigid support to the lumbar spine.
- Stabilized the back and enhanced faster recovery.
- Adjustable straps to provide maximum support and comfortability.

Feature

Application

Low back pain
Lumbar spondylosis
Spondylolisthesis
Back strain

Material

Polyester, Cotton, Rubber, Nylon,
Steel stay.

Quantity 1 PC

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in
XXL	109 - 119	cm	43 - 47	in





CPO-6204

Back Support Auxi

Breathable lumbar support

- Stabilized lumbar spine
- 8 inches height with four steel stays to provide maximum support to the back.
- Slim breathable mesh decrease muggy feeling.

Feature

Application

Low back pain
Lumbar spondylosis
Back strain


Material

Nylon, Polyester, Cotton, Rubber,
Steel stay.

Quantity 1 PC

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in





CPO-6206

Back Support ST

Magnetic back belt

- Magnet design in the waist and belly help to increase circulation.
- Flexible and keep mobility unrestricted. Suitable for health care use.
- 6 inches height to support lower back.

Feature

Application

Low back pain
Back strain

Material

Polyester, Nylon, Cotton, Rubber,
Magnet.

Quantity 1 PC

S	58 - 68	cm	23 - 27	in
M	68 - 80	cm	27 - 31 ½	in
L	80 - 92	cm	31 ½ - 36	in
XL	92 - 105	cm	36 - 41 ¼	in





CPO-6207

Back Support ST+

Lumbar belt

- 6 inches height, lightweight with an elastic strap that makes it fit perfectly.
- Two plastic stays provide mild support to the lower back.
- Increase the abdominal pressure to promote muscle strength in the waist.

Feature

Application

Low back pain
Back strain

Material

Cotton, Nylon, Rubber.

Quantity 1 PC

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



Back Support TL

6" high lumbar support



CPO-6208

Feature

- Comfortable and lightweight to carry around the waist.
- Comfortable to wear during the activities.
- Adjustable straps to increase perfect fit and compression to the lower back.
- Four plastic stays provide slightly support in waist.

Material

Cotton, Nylon, Polyester, Rubber, PP.

Quantity

1 PC

Application

Low back pain
Lumbar spondylosis
Back strain

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



Back Support Auxi+

Breathable lumbar support



CPO-6211

Feature

- 9 inches height with four steel stays provides solid protection and support of the lower back.
- Slim breathable mesh decreases muggy feeling.
- Adjustable strengthen straps increase perfectly fit and compression of the belt.

Material

Nylon, Polyester, Rubber, Steel stay.

Quantity

1 PC

Application

Low back pain
Lumbar spondylosis
Back strain

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



Lumbar Support with Suspenders

Industrial back support



CPO-6212

Feature

- 8 inches height with four plastic stays provide stable and flexible support of back.
- Adjustable strengthen straps increase perfect fit and compression to the lumbar spine.
- Adjustable shoulder suspension straps.

Material

Nylon, Polyester, Rubber, Steel stay.

Quantity

1 PC

Application

Low back pain
Back strain

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



Orthopedic Support



Back Support AG

Breathable lumbar support

Feature

- 8 inches height with two plastic stays provide stable and flexible support of back.
- Slim breathable mesh material decrease muggy feeling.
- Comfortable to wear during the activities.
- Adjustable strengthen straps increase the perfectly fit and compression of the belt.

Material

Nylon, Polyester, Rubber, Cotton, PP.

Quantity

1 PC

Application

Low back pain
Back strain

S	58 - 68	cm	23 - 27	in
M	68 - 80	cm	27 - 31 ½	in
L	80 - 92	cm	31 ½ - 36	in
XL	92 - 105	cm	36 - 41 ¼	in



Back Support GT

Lumbar support

Feature

- 9 inches height with four steel stays to provide solid support to the lumbar spine.
- Adjustable straps to increase perfect fit and suitable compression.
- Instant pain relief and enhanced faster recovery.

Material

Nylon, Polyester, Rubber, Steel stay.

Quantity

1 PC

Application

Low back pain
Lumbar spondylosis
Back strain
Spondylolisthesis

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



Back Support UE

Breathable lumbar support

Feature

- 9 inches height with four plastic stays provides slightly support in waist.
- Slim breathable mesh material decrease muggy feeling.
- Comfortable to wear during the activities.

Material

Nylon, Polyester, Rubber, Cotton, PP, PVC, PE.

Quantity

1 PC

Application

Low back pain
Lumbar spondylosis
Back strain

S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in
XXL	109 - 119	cm	43 - 47	in





AerGo Elbow Support

Elbow support

Feature

- Protect the weak elbow.
- Provide proper amount of pressure to the joint.
- 3D airknit technology.

Application

Weak elbow joint
Elbow pain

Material

Spandex, Nylon.

Quantity

1 PC



S	15 - 20	cm	6 - 8	in
M	20 - 25	cm	8 - 10	in
L	25 - 30	cm	10 - 12	in
XL	30 - 35	cm	12 - 14	in

CPO-7305



AerGo Knee Support TL

Knee support

Feature

- Protect the weak knee.
- Provide proper amount of pressure to the joint.
- 3D airknit technology.

Application

Weak knee joint
Knee pain

Material

Spandex, Nylon.

Quantity

1 PC



S	28.5 - 33.5	cm	11 ¼ - 13 ¼	in
M	33.5 - 39	cm	13 ¼ - 15 ½	in
L	39 - 44.5	cm	15 ½ - 17 ½	in
XL	44.5 - 50	cm	17 ½ - 19 ¾	in

CPO-7601



AerGo Knee Support TL+

Knee support with spiral stay

Feature

- Spiral stays added to stabilized weak knee.
- Provide proper amount of pressure to the joint.
- 3D airknit technology.

Application

Osteoarthritis
Weak knee joint
Knee pain

Material

Spandex, Nylon, Steel.

Quantity

1 PC



S	28.5 - 33.5	cm	11 ¼ - 13 ¼	in
M	33.5 - 39	cm	13 ¼ - 15 ½	in
L	39 - 44.5	cm	15 ½ - 17 ½	in
XL	44.5 - 50	cm	17 ½ - 19 ¾	in

CPO-7602



AerGo Ankle Support

Ankle support

Feature

- Protect the weak ankle.
- Provide proper amount of pressure to the joint.
- 3D airknit technology.

Application

Weak ankle
Ankle pain

Material

Spandex, Nylon.

Quantity

1 PC



S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in

CPO-7702



CPO-7202

AerGo Back Brace

Lumbar support

Feature

- 9 inches height.
- 4 steel stays provide stable support of back.
- Breathable mesh fabric decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

Material

Polyester, Nylon, Rubber, PP, Steel.

Application

Low back pain
Lumbar spondylosis
Back strain

Quantity

1 PC



S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



CPO-7203

Reinforced Back Support

Reinforced Back Support

Feature

- 9 Inches height.
- 4 Steel stays provide stable support of the back.
- Breathable mesh fabric decreases the muggy feeling.
- Outside strengthen strap whose design is easy for users to pull increases perfectly fit and compression of the belt.

Material

Polyester, Nylon, Rubber, PP, Steel, PVC.

Application

Low back pain
Lumbar spondylosis
Back strain

Quantity

1 PC



S	58 - 70	cm	23 - 27 ½	in
M	70 - 85	cm	27 ½ - 33 ½	in
L	85 - 100	cm	33 ½ - 39 ½	in
XL	100 - 110	cm	39 ½ - 43 ¼	in



CPO-7704

AerGo Toe Splint for Hallux Valgus

Toe splint for Hallux valgus

Patented

Feature

- Bilateral light aluminum hinges serve as the prop of toes and keep them fixed.
- Detachable hinge on little toes can be applied if need.
- Relieves discomfort caused by hallux valgus and reduces friction between toes.
- Breathable mesh fabric for comfort feeling.
- Can be used while walking.

Material

Neoprene, Nylon, Aluminum, Polyester.

Application

Hallux valgus
Tailor's bunion
Varus toe

Quantity

1 PC



S	20 - 24	cm	8 - 9 ½	in
L	24 - 29	cm	9 ½ - 11 ½	in

Aergo Golf / Tennis Elbow Support

Tennis / Golf elbow support

Patented

Feature

- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.
- Silicone pad provides extra compression to stabilize the injured area.
- Light weight design helps to increase the performance for sport.
- Easy to wear and adjustable.

Application

Tennis elbow
Golfer's elbow

Material

Neoprene, Nylon, Silicone pad,
Polyester.

Quantity

1 PC

S	20 - 28	cm	8 - 11	in
L	28 - 34	cm	11 - 13 ½	in



CPO-7304

Aergo Wrist Brace

Wrist brace

Feature

- For carpal tunnel syndrome use.
- Aluminum stay enhance the stability of wrist joint.
- Easy to wear and adjust.

Application

Carpal tunnel syndrome
Wrist sprain
Colles fracture

Material

Polyester, Nylon, Aluminum stay.

Quantity

1 PC

S	14.5 - 16.5	cm	5 ¾ - 6 ½	in
M	16.5 - 18.5	cm	6 ½ - 7 ¼	in
L	18.5 - 21.5	cm	7 ¼ - 8 ½	in
XL	21.5 - 23.5	cm	8 ½ - 9 ¼	in



CPO-7401

Aergo Ankle Brace

Ankle brace

Feature

- Assist the rehabilitation of ankle sprain.
- Two plastic pads enhance the stability.
- For ankle joint stabilization.

Application

Ankle instability
Ankle sprain / strain

Material

Nylon, Polyester, PE.

Quantity

1 PC

S	15.5 - 21	cm	6 - 8 ¼	in
M	21 - 26	cm	8 ¼ - 10 ¼	in
L	26 - 31	cm	10 ¼ - 12 ¼	in
XL	31 - 36	cm	12 ¼ - 14 ¼	in



CPO-7701

Aergo Knee Support

Elastic knee support

Feature

- Bilateral spiral stays provide stable support for knee joint.
- Breathable elastic material decrease muggy feeling.
- To maintain the structure of knee joint and provide suitable compression on the knee joint.

Application

Patella instability
Knee weakness
Knee pain

Material

Rayon, Polyester, Steel stay.

Quantity

1 PC

S	27 - 30	cm	10 ½ - 11 ¾	in
M	30 - 34	cm	11 ¾ - 13 ¼	in
L	34 - 40	cm	13 ¼ - 15 ¾	in
XL	40 - 44	cm	15 ¾ - 17 ¼	in
XXL	44 - 48	cm	17 ¼ - 19	in



CPO-6601



Graduated compression design helps improve better blood flow and circulation.



Integral knitting design for heel and toe provides comfortable and easy wearing.



Patent Lycra material provides soft texture, suitable for long time wearing.



Elastic band prevents slippery and binding.



Open toe design helps monitor blood circulation.



Soft band with silicone dots prevents slippery and binding.



Special design with less elasticity helps shape up hips and improves the circulation of legs.

Medical Compression Stockings series

Medical compression stockings help fight blood pressure within the veins, mitigate or avoid aggravation of damage. Graduated compression design helps improve blood flow and circulation, and it can help the blood return and increase the effects of circulation metabolism.

Features of Compression Stockings

- Design for German regulatory compression standard, ccl.1:18-21mmHg, ccl.2:23-32mmHg, and ccl.3:34-46mmHg. The graduated compression design helps improve blood flow and circulation, reduce the edema and aching of lower extremities.
- Graduated compression keeps legs in good conditions, and effectively eliminates the edema and shapes up the legs.
- Lycra yarn has better elasticity and durability, and it can maintain product efficacy and stand the wear and tear. The soft texture provides users with long-time wearing comforts.
- Non-toxic yarn reduces allergy symptoms. The yarn is granted certificate of Oeko-Tex Standard 100 which reassures customers of its premium quality.
- Ccl.1 compression stockings are suitable for people need to endure long-term standing such as for teachers, sales assistants, waiters, flight attendants, and people who have edema or fatigued, aching lower limbs. Those who suffer from slight varicosity can use ccl.1 to prevent varicosity from getting worse. People with moderate varicosity or post-surgery patients can use ccl.3 for restoration, and use ccl.2 for the long term maintenance.

Applications of each compression level

ccl.1:18-21mmHg

- Minor varicosity.
- Tired, aching legs.
- Minor ankle, leg and foot swelling.
- Reduce leg swelling and prevent varicosity during pregnancy.

ccl.2:23-32mmHg

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

ccl.3:34-46mmHg

- Severe varicosity.
- Severe edema.
- Post surgical use.

Resile Series

CCLI:18-21mmHg | CCLII:23-32mmHg | CCLIII:34-46mmHg



Resile™ Urban
Knee High closed toe

Application

- Mild to severe varicosity or spider veins.
- Mild to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Prevent varicosity and leg swelling during pregnancy.
- Post-surgical use.

Material
75% Nylon, 25% Lycra.

Quantity 1 PR



Size	1 - 7
Color	 
Suitable for	 

Class I	CPS-3001
Class II	CPS-3002
Class III	CPS-3006



Resile™ Soft
Knee High open toe

Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material
75% Nylon, 25% Lycra.

Quantity 1 PR



Size	1 - 7
Color	 
Suitable for	 

Class II	CPS-3004
Class III	CPS-3007



Resile™ Fit
Mid Thigh closed toe

Application

- Mild to severe varicosity or spider veins.
- Mild to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Prevent varicosity and leg swelling.
- Post-surgical use.

Material
75% Nylon, 25% Lycra.

Quantity 1 PR



Size	1 - 7
Color	 
Suitable for	 

Class I	CPS-3101
Class II	CPS-3102
Class III	CPS-3105



Resile™ Sophisticate
Mid Thigh open toe

Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material
75% Nylon, 25% Lycra.

Quantity 1 PR



Size	1 - 7
Color	 
Suitable for	 

Class II	CPS-3104
Class III	CPS-3106

CCLI:18-21mmHg | CCLII:23-32mmHg | CCLIII:34-46mmHg

Resile™ Dynamic

Thigh High closed toe



Application

- Minor to severe varicosity or spider veins.
- Minor to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Post-surgical use.

Material

75% Nylon, 25% Lycra.

Quantity

1 PR



Size	1 - 7
Color	
Suitable for	

Class I	CPS-3301
Class II	CPS-3302
Class III	CPS-3305

Resile™ Flexi

Thigh High open toe



Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

75% Nylon, 25% Lycra.

Quantity

1 PR




Size	1 - 7
Color	
Suitable for	

Class II	CPS-3304
Class III	CPS-3306

Resile™ Smooth

Pantyhose closed toe



Application

- Minor to severe varicosity or spider veins.
- Minor to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Post-surgical use.

Material

75% Nylon, 25% Lycra.

Quantity

1 PR




Size	1 - 7
Color	
Suitable for	

Class I	CPS-3401
Class II	CPS-3402
Class III	CPS-3405

Resile™ Comfy

Pantyhose open toe



Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

75% Nylon, 25% Lycra.

Quantity

1 PR



Size	1 - 7
Color	
Suitable for	

Class II	CPS-3404
Class III	CPS-3406

Euniiz Series

CCLI:18-21mmHg | CCLII:23-32mmHg | CCLIII:34-46mmHg



Euniiz™ Urban

Knee High closed toe

Application

- Mild to severe varicosity or spider veins.
- Mild to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Prevent varicosity and leg swelling during pregnancy.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR



Size	1-5	1-7
Color		
Suitable for		

Class I	CPS-2001
Class II	CPS-2002
Class III	CPS-2006



Euniiz™ Soft

Knee High open toe

Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR



Size	1	-	7
Color			
Suitable for			

Class II	CPS-2004
Class III	CPS-2007



Euniiz™ Fit

Mid Thigh closed toe

Application

- Mild to severe varicosity or spider veins.
- Mild to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Prevent varicosity and leg swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR



Size	1-5	1-7
Color		
Suitable for		

Class I	CPS-2101
Class II	CPS-2102
Class III	CPS-2105



Euniiz™ Sophisticate

Mid Thigh open toe

Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR



Size	1	-	7
Color			
Suitable for			

Class II	CPS-2104
Class III	CPS-2106

CCLI:18-21mmHg | CCLII:23-32mmHg | CCLIII:34-46mmHg

Euniiz™ Dynamic

Thigh High closed toe



Application

- Minor to severe varicosity or spider veins.
- Minor to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR



Size	1-5	1-7
Color		
Suitable for	♀	♂

Class I	CPS-2301
Class II	CPS-2302
Class III	CPS-2201

Euniiz™ Flexi

Thigh High open toe



Application

- Moderate to severe varicosity or spider veins
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR




Size	1 - 7	
Color		
Suitable for		

Class II	CPS-2304
Class III	CPS-2202

Euniiz™ Smooth

Pantyhose closed toe



Application

- Minor to severe varicosity or spider veins.
- Minor to severe edema.
- Tired, aching legs.
- Reduce ankle, leg, and foot swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity

1 PR




Size	1-5	1-7
Color		
Suitable for	♀	

Class I	CPS-2401
Class II	CPS-2402
Class III	CPS-2405

Euniiz™ Comfy

Pantyhose open toe



Application

- Moderate to severe varicosity or spider veins.
- Moderate to severe edema.
- Reduce leg swelling.
- Post-surgical use.

Material

80% Nylon, 20% Lycra.

Quantity


1 PR



Size	1 - 7	
Color		
Suitable for		

Class II	CPS-2404
Class III	CPS-2406

Medical Compression Stockings



Size		Compression stockings						
------	--	-----------------------	--	--	--	--	--	--



Size		1	2	3	4	5	6	7
Circumference	cH	≤ 95 cm ≤ 37½"	≤ 100 cm ≤ 39¼"	≤ 105 cm ≤ 41¼"	≤ 110 cm ≤ 43¼"	≤ 114 cm ≤ 45"	≤ 118 cm ≤ 46½"	≤ 120 cm ≤ 47¼"
	cG	43 - 48 cm 17" - 18¾"	45 - 52 cm 17¾" - 20½"	47 - 56 cm 18½" - 22"	50 - 60 cm 19¾" - 23¾"	53 - 63 cm 20¾" - 24¾"	56 - 66 cm 22" - 26"	58 - 68 cm 22¾" - 26¾"
	cC	28 - 34 cm 11" - 13¼"	30 - 37 cm 11¾" - 14½"	33 - 40 cm 13" - 15¾"	35 - 43 cm 13¾" - 17"	37 - 45 cm 14½" - 17¾"	39 - 47 cm 15¼" - 18½"	41 - 49 cm 16" - 19¼"
	cB	18 - 20 cm 7" - 7¾"	20 - 22 cm 7¾" - 8¾"	22 - 24 cm 8¾" - 9½"	24 - 26 cm 9½" - 10¼"	26 - 28 cm 10¼" - 11"	28 - 30 cm 11" - 11¾"	30 - 32 cm 11¾" - 12½"
Size		1	2	3	4	5	6	7
Length	IG	62 - 70 cm 24½" - 27½"	62 - 70 cm 24½" - 27½"	66 - 74 cm 26" - 29"	66 - 74 cm 26" - 29"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"
	ID	30 - 38 cm 11¾" - 15"	30 - 38 cm 11¾" - 15"	34 - 40 cm 13¼" - 15¾"	34 - 40 cm 13¼" - 15¾"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"

cH The circumference of hip. **cG** The circumference of thigh. **cC** The circumference of calf.
cB The circumference of ankle. **IG** The length from heel to hip. **ID** The length from heel to knee.

Anti-embolism stockings

Knee high anti-embolism stockings 18mmHg



Application

- For post-surgical, reclining or bed-confined patients to get better circulation.
- Reduce the risk of embolism formation.
- Prevent thrombo-embolic disease.

Material

Nylon, Spandex.

Quantity 1 PR

Size 1 - 7

Color

CPS-2005



Thigh high anti-embolism stockings with silicone top band 18mmHg



Application

- For post-surgical, reclining or bed-confined patients to get better circulation.
- Reduce the risk of embolism formation.
- Prevent thrombo-embolic disease.

Material

Nylon, Spandex.

Quantity 1 PR


Size 1 - 7

Color

CPS-2309



Size | Anti-embolism stockings



Size		1	2	3	4	5	6	7
Circumference	cG	43 - 48 cm 17" - 18 3/4"	45 - 50 cm 17 3/4" - 19 3/4"	47 - 52 cm 18 1/2" - 20 1/2"	50 - 56 cm 19 3/4" - 22"	53 - 58 cm 20 3/4" - 22 3/4"	56 - 60 cm 22" - 23 1/2"	58 - 64 cm 22 3/4" - 25 1/4"
	cC	28 - 34 cm 11" - 13 1/4"	30 - 37 cm 11 3/4" - 14 1/2"	33 - 40 cm 13" - 15 3/4"	35 - 43 cm 13 3/4" - 17"	37 - 45 cm 14 1/2" - 17 3/4"	39 - 47 cm 15 1/4" - 18 1/2"	41 - 49 cm 16" - 19 1/4"
	cB	18 - 20 cm 7" - 7 3/4"	20 - 22 cm 7 3/4" - 8 3/4"	22 - 24 cm 8 3/4" - 9 1/2"	24 - 26 cm 9 1/2" - 10 1/4"	26 - 28 cm 10 1/4" - 11"	28 - 30 cm 11" - 11 3/4"	30 - 32 cm 11 3/4" - 12 1/2"
Size		1	2	3	4	5	6	7
Length	IG	62 - 70 cm 24 1/2" - 27 1/2"	62 - 70 cm 24 1/2" - 27 1/2"	66 - 74 cm 26" - 29"	66 - 74 cm 26" - 29"	68 - 76 cm 26 3/4" - 30"	68 - 76 cm 26 3/4" - 30"	68 - 76 cm 26 3/4" - 30"
	ID	30 - 38 cm 11 3/4" - 15"	30 - 38 cm 11 3/4" - 15"	34 - 40 cm 13 1/4" - 15 3/4"	34 - 40 cm 13 1/4" - 15 3/4"	36 - 42 cm 14 1/4" - 16 1/2"	36 - 42 cm 14 1/4" - 16 1/2"	36 - 42 cm 14 1/4" - 16 1/2"

cG The circumference of thigh. cC The circumference of calf. cB The circumference of ankle.
IG The length from heel to hip. ID The length from heel to knee.

Medical Compression Stockings

Travel Stockings



Knee high travel stockings 18-21 mmHg

Application

- Minor ankle, leg and foot swelling during travel.
- Helps prevent DVT.

Material

Nylon, Spandex.


Quantity

1 PR



Size	S	-	XL
Color		-	
Suitable for		-	
CPS-1002			

Size | Travel stockings




Size		S	M	L	XL
Circumference	cC	28 - 34 cm 11" - 13 ¼"	33 - 40 cm 13" - 15 ¾"	37 - 45 cm 14½" - 17¾"	41 - 49 cm 16" - 19 ¼"
	cB	18 - 20 cm 7" - 7¾"	22 - 24 cm 8¾" - 9½"	26 - 28 cm 10¼" - 11"	30 - 32 cm 11¾" - 12½"
Size		S	M	L	XL
Length	ID	30 - 38 cm 11¾" - 15"	34 - 40 cm 13¼" - 15¾"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"

cC The circumference of calf. cB The circumference of ankle. ID The length from heel to knee.

Diabetic Socks



Diabetic Socks



Feature

- Antibacterial silver fiber to prevent the growing of bacteria.
- Seamless stitching to prevent the friction on toes.
- Elastic and pressure free top band to avoid oppression.
- Comfortable toweling sole provides cushion and protection for feet.

Application

- Health caring for diabetic.
- Poor circulation in feet.
- Smelly feet.

Material

Cotton, Nylon, Silver fiber, Lycra.

Quantity

1 PR

Size	S - L	
Color	☒	☐
Suitable for	☒	☐
CPS-4001		

SIZE	EUR	UK	US	
			MEN	WOMEN
S	34 - 37	1.5 - 4	3½ - 5½	4½ - 6½
M	38 - 42	4 - 7.5	4½ - 8	6 - 9½
L	43 - 47	6 - 10	6½ - 10½	8 - 11½



Insoles and Footcare

Keep Comfort Health Plus

- Medical grade silicone minimizes allergy reactions.
- Medical grade silicone never flattens, deforms or degrades for long time use.
- Gently massage the feet.
- Helps relieve foot discomfort and fatigue from prolonged standing or walking.
- Redistributes pressure on the foot.
- The ergonomic design helps wearing with ease.
- Silicone and gel materials are washable and reusable



Medical grade silicone minimizes allergy reactions.



Multiple density design helps release and redistribute pressure on the feet.



The ergonomic design helps wearing with ease.



Prevent uncomfortable friction.



Skin moisturizing.



Slim and lightweight design.

Insoles and Footcare

Silicone

Silicone Insoles



- Feature
- Softer inserts provides more cushioning and shock absorption on metatarsal and heel areas.
 - Metatarsal pad helps to support metatarsal arch.
 - Multiple holes for air circulation.


- Application
- Heel pain.
Metatarsalgia.
Relieves the discomfort of plantar fasciitis.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
34	33 - 34	1 - 2	2 - 3	3½ - 4½
36	35 - 36	2½ - 3½	3½ - 4½	5 - 6
38	37 - 38	4 - 5	5 - 6	6½ - 7½
40	39 - 40	5½ - 6½	6½ - 7½	8 - 9
42	41 - 42	7 - 8	8 - 9	9½ - 10
44	43 - 44	8½ - 9½	9½ - 10½	10½ - 11½
46	45 - 46	10 - 11	11 - 12	12 - 13½
48	47 - 48	11½ - 13	12½ - 13½	14 - 15

Silicone insoles with toe support



- Feature
- Softer inserts provides more cushioning and shock absorption on lateral of mid foot and heel areas.
 - Raised toe support helps to support toe area and relieve the discomfort of toes.


- Application
- Heel pain.
Relieves the pressure in lateral side.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	36 - 38	3½ - 5	4½ - 6	6 - 7½
M	39 - 41	5½ - 7½	6½ - 8½	8 - 9½
L	42 - 44	8 - 9½	9 - 10½	10 - 11½

Silicone insoles



- Feature
- Soft textile liner to provide extra comfort.
 - Can be cut to fit the shoe size.


- Application
- Foot pain.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	34 - 37	2 - 4½	3 - 5½	4½ - 7
M	38 - 42	5 - 8	6 - 9	7½ - 10
L	43 - 47	8½ - 11½	9½ - 12½	10½ - 14

Silicone

Silicone insoles



Feature

- Softer inserts provides more cushioning and shock absorption on metatarsal and heel areas.
- Metatarsal pad helps to support metatarsal arch.

Application


Heel pain.
Relieves the discomfort of
plantar fasciitis.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
XS	37 - 38	4 - 5	5 - 6	6½ - 7½
S	39 - 40	5½ - 6½	6½ - 7½	8 - 9
M	41 - 42	7 - 8	8 - 9	9½ - 10
L	43 - 44	8½ - 9½	9½ - 10½	10½ - 11½
XL	45 - 46	10 - 11	11 - 12	12 - 13½

Silicone insoles



Feature

- Softer inserts provides more cushioning and shock absorption on arch and heel area.
- Slim and light weight design.

Application


Heel pain.
Relieves the pressure in medial
side.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 40	4 - 6½	5 - 7½	6½ - 9
L	41 - 44	7 - 9½	8 - 10½	9½ - 11½

Silicone insoles ¾ length



Feature

- Softer inserts provides more cushioning and shock absorption on metatarsal and heel areas.
- Provides heel stability and supports the longitudinal arch.

Application


Heel pain.
Metatarsalgia.
Relieves the discomfort of
plantar fasciitis.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 40	4 - 6½	5 - 7½	6½ - 9
L	41 - 44	7 - 9½	8 - 10½	9½ - 11½

Insoles and Footcare

Silicone

Silicone heel cups



Feature

- Softer insert provides more cushioning and shock absorption on heel area.
- Cup design provides heel stability and comfort.

Application


Heel Pain.
Applicable for long time standing.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11
XL	44 - 46	9 - 11	10 - 12	11½ - 13½

Silicone heel pads with removable pads



Feature

- Helps absorb impact and prevent heel pain.
- Removable pads helps relieve heel spur.
- Could be used with/without removable pads for different needs.

Application


Heel Pain.
Heel spur.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11

Silicone heel pads



Feature

- Softer insert provides more cushioning and shock absorption on heel area.
- Redistributes pressure on the foot.
- Non-slip design on the heel area.

Application


Heel Pain.
Applicable for long time standing.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11

Silicone

Silicone heel cushions



CPF-1204

Feature

- Softer insert on heel area helps absorb shock and relieve pain.
- Redistributes pressure on the foot.
- Different thickness of medial and lateral edge could be used for pronation or supination.

Application


Heel Pain.
Applicable for long time standing.
Slight varus or valgus.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	≤ 34	≤ 2	≤ 3	≤ 4½
M	35 - 40	2½ - 6½	3½ - 7½	5 - 9
L	41 - 46	7 - 11	8 - 12	9½ - 13½

Silicone heel pads



CPF-1205

Feature

- Squares design helps absorb shock and relieve pain on heel area.
- Redistributes pressure on the foot.
- Thin type fits most shoe types.

Application


Heel Pain.
Applicable for long time standing.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11

Silicone heel spur pads



CPF-1206

Feature

- U-shaped design helps shock absorption and pressure relief on the heel area.
- Provides heel stability and prevents ankle sprain.
- Relieves heel pain and heel spur.
- Easy to fit the shoes.

Application


Heel Pain.
Calcaneus instability.
Flatfoot.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11

Insoles and Footcare

Silicone

Silicone heel cups



CPF-1208

Feature

- Soft material provides cushioning and shock absorption on heel area.
- Cup design provides heel stability and prevents ankle sprain.
- Ergonomic design fits the shape of foot, provides comfortable wearing.

Application

Heel Pain.
Applicable for long time standing.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11

Silicone metatarsal pads



CPF-1301

Feature

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.
- Helps realign toes.
- Toe spreader design could be used for hallux valgus.

Application


Metatarsalgia.
Forefoot calluses.
Overlapping toes.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½

Silicone metatarsal pads



CPF-1302

Feature

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.

Application

Metatarsalgia.
Forefoot calluses.

Material

Silicone.

Quantity

1 PR



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½

Silicone

Silicone toe caps



CPF-1401

Feature

- Provides cushioning and protection to toes.
- Relieves discomfort caused by calluses and corns.
- Fit your toes or fingers.

Application


Overlapping toes.
Calluses on toe.
Corns on toe.

Material

Silicone.

Quantity

2 PCS



S

Diameter 16mm

L

Diameter 20mm



Silicone toe tubes



CPF-1402

Feature

- Provides cushioning and protection to toes.
- Relieves discomfort caused by calluses and corns.
- Separate and straighten overlapping toes.

Application


Overlapping toes.
Calluses on toe.
Corns on toe.

Material

Silicone.

Quantity

2 PCS



S

Diameter 15mm

M

Diameter 18mm

L

Diameter 22mm



Silicone toe crests



CPF-1501

Feature

- Helps reduce pressure and friction.
- Stretch contracted toes.
- Relieves discomfort caused by hammer, claw, and Mallet toes.
- Toe loop design keeps toe crest in place.

Application

Hammer toes.
Claw toes
Mallet toes

Material

Silicone.

Quantity

1 PR



SIZE

EUR

UK

US

MEN

WOMEN

S

33 - 36

1 - 3½

2 - 4½

3½ - 6

M

37 - 41

4 - 7

5 - 8

6½ - 9½

L

42 - 46

7½ - 11

8½ - 12

10 - 13½



Insoles and Footcare

Silicone

Silicone toe spreaders



CPF-1502

Feature

- Helps straighten and realign toes.
- Relieves discomfort caused by overlapping toes and hallux valgus.
- Reduces friction between toes.
- Separates the overlapping toes.

Application


Hallux valgus.
Overlapping toes.

Material

Silicone.

Quantity

2 PCS



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

Silicone toe separators



CPF-1503

Feature

- Helps straighten and realign toes.
- Relieves discomfort caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

Application


Hallux valgus.
overlapping toes.

Material

Silicone.

Quantity

2 PCS



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

Silicone tailor's bunion protector



CPF-1504

Feature

- Helps reduce pressure and friction on 5th toe area.
- Toe loop design keeps protector in place.

Application

Tailor's bunion.
Calluses.

Material

Silicone.

Quantity

1 PC



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½

Silicone

Silicone hallux valgus protector




- Feature**
- Helps reduce pressure and friction on hallux area.
 - Toe loop design keeps protector in place.

Application Hallux valgus.
Bunion.

Material Silicone

Quantity 1 PC



SIZE	EUR	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½

Silicone corn pad



- Feature**
- Soft silicone pad protects and prevents of corns and callous formation.

Application Calluses.
Corns.

Material Silicone, Polyester.

Quantity 4 PCS



One Size

Silicone corn pads




- Feature**
- Soft silicone pad protects and prevents of corns and callous formation.

Application Calluses.
Corns.

Material Silicone, Polyester.

Quantity 9 PCS



One Size

Silicone heel liners




- Feature**
- Prevents the friction between shoes and back of the heel.

Application Back of heel frictions.
Calluses.

Material Silicone

Quantity 1 PR



One Size

Insoles and Footcare

High heel & Sandals

High heel silicone insoles



- Metatarsal pad helps to support metatarsal arch.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

Application Suitable for high heel shoes. Metatarsalgia.

Material Silicone

Quantity 1 PR



SIZE	EUR	UK	US
S	35 - 37	2½ - 4	5 - 6½
M	38 - 40	4½ - 6½	7 - 9
L	41 - 43	7 - 8½	9½ - 11

High heel silicone insoles



- Softer inserts provide more cushioning and shock absorption on metatarsal area.
- Metatarsal pad helps to support metatarsal arch.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

Application Suitable for high heel shoes. Metatarsalgia.

Material Silicone.

Quantity 1 PR



SIZE	EUR	UK	US
S	35 - 37	2½ - 4	5 - 6½
M	38 - 40	4½ - 6½	7 - 9
L	41 - 43	7 - 8½	9½ - 11

High heel Silicone metatarsal pads



- Metatarsal pad helps to support metatarsal arch.
- Slim and light weight fits most women's fashion footwear and high heel shoes style.

Application Suitable for high heel shoes. Metatarsalgia.

Material Silicone.

Quantity 1 PR



SIZE	EUR	UK	US
S	33 - 37	1 - 4	3½ - 6½
L	38 - 42	5 - 8	7 - 10

High heel Silicone metatarsal pads with removable dot



- Metatarsal pad helps to support metatarsal arch.
- Removable dot could be adjusted for different support region.
- Slim and light weight fits most women's fashion footwear and high heel shoes style.

Application Suitable for high heel shoes. Metatarsalgia.

Material Silicone

Quantity 1 PR



SIZE	EUR	UK	US
S	33 - 37	1 - 4	3½ - 6½
L	38 - 42	5 - 8	7 - 10

Other

Gel hand therapy ball (Round) CPF-2604

Gel hand therapy ball (Egg) CPF-2605



Application Hand rehabilitation use.
Hand exercise use.

Material Gel.

Quantity 1 PC



Pink	Extra soft
Purple	Soft
Green	Medium
Orange	Hard
Black	Extra hard

Gel

Gel bunion protector with toe spreader




Feature

- Helps straighten and realign toes.
- Relieves discomfort caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

Application Hallux valgus.
Bunion.
Overlapping toes.

Material Gel.

Quantity 1 PC



One Size



Insoles and Footcare

Moisturizing

Moisturizing socks



CPF-8601

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
 - Moisturizing gel improves dry skin.
 - Lavender essential oil helps recovering weary spirit.
- Calluses and stratum corneum.
Dry skin.

Feature

Application

Material

Outer sleeve: Nylon, Spandex.
Inner gel: Gel, Jojoba oil, Lavender essential oil.

Quantity 1 PR



One Size

Moisturizing heel socks



CPF-8602

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
 - Moisturizing gel improves dry skin.
 - Lavender essential oil helps recovering weary spirit.
- Calluses and stratum corneum.
Dry skin.

Feature

Application

Material

Outer sleeve: Nylon, Spandex.
Inner gel: Gel, Jojoba oil, Lavender essential oil.

Quantity 1 PR



One Size

Moisturizing gloves



CPF-8603

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
 - Moisturizing gel improves dry skin.
 - Lavender essential oil helps recovering weary spirit.
- Calluses and stratum corneum.
Dry skin.

Feature

Application

Material

Outer sleeve: Nylon, Spandex.
Inner gel: Gel, Jojoba oil, Lavender essential oil.

Quantity 1 PR



One Size

Moisturizing elbow sleeves



CPF-8604

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
 - Moisturizing gel improves dry skin.
 - Lavender essential oil helps recovering weary spirit.
- Calluses and stratum corneum.
Dry skin.

Feature

Application

Material

Outer sleeve: Nylon, Spandex.
Inner gel: Gel, Jojoba oil, Lavender essential oil.

Quantity 1 PR



One Size

COMFORT PLUS CORPORATION

ADD | 11F., NO 123-9, Xingde Rd., Sanchong District, New Taipei City 24158, Taiwan, R.O.C.

TEL : +886-2-85113758 FAX : +886-2-85113755

E-mail | service@comfortplus-medical.com

<http://www.comfortplus-medical.com>

<https://www.eunicemed.com>